



#sharing challenges
and solutions in practice

FUJIFILM Diosynth Biotechnologies CDMO: Multifunctional Single-Use Purification System for connected and integrated Continuous Processing

Facility & Technology Projects

Part of PharmaCongress – Düsseldorf/Neuss, 31 May - 1 June 2022

SYMPHON X™

Single Use Technology collaboration case study

Development and use of a multi-
functional purification skid

FUJIFILM Diosynth Biotechnologies
Gemü GmbH




FUJIFILM Diosynth Biotechnologies



www.fujifilmdiosynth.com

End-to-End CDMO



www.fujifilmdiosynth.com

Strategies for Supplying Market Demand



Scale-Up

Multiple 20,000-L
bioreactors

Efficiency of scale

Ideal for high-volume
products



Scale-Out

Multiple 2,000 L bioreactors
Options for multiplexing (2 x
2000L)

Flexible strategy mitigates
risks due to uncertainties in
commercial demand

Ideal for lower volume
products

Leverages FDB's 'mAb
platform' technology



Continuous

Single or multiple 500-L
bioreactors or larger

Improved production
efficiencies


Flexibility and potential
to deliver high-volume

www.fujifilmdiosynth.com

Made by GEMÜ

Our quality promise

- Manufacturing in accordance with internationally recognised **production and security standards**
- **A consistently high level of quality all over the world** – irrespective of the manufacturing site



Why SymphonX™?


Normally in the DSP several units are needed for different processing steps:

- Tangential flow filtration
- Protein A affinity chromatography
- Viral Inactivation
- Cation and anion exchange chromatography
- Viral Filtration
- Ultrafiltration/diafiltration

This requires a lot of equipment types and various training needs!



Successful Collaboration



Description

SymphonX™

- Purification platform system
- Integrated buffer dilution system
- Single flow path for multiple use
- Versatile modular system
- Scalable system 6 – 740 L/hr



Key Features

SymphonX™

- Single-use flow path with disposable pump head, valve bodies and sensor probes
- Flexible software interface
- Automation friendly process control
- Developed and validated for use in GMP
- Successfully deployed for
 - GMP batch manufacturing
 - non-GMP continuous mAb production



Main Advantages

SymphonX™

- Multiple application capabilities within the same skid i. e. one skid can cover the whole Bio purification train
- Automated system
- Ease of training (one skid as opposed to many)
- Tech transfer between locations and facilities
- Improved forecastable supply chain
- Embedding sustainability principles into Single-Use Technology through robust flow path design




User Requirement Specification

SymphonX™

- Sustaining a linear flow within a range of 6 – 740 L/hr
- Internal diameter of all flow path components is ½”
- Materials PP, EPDM, PEEK, Silicon, TPE or PSU
- Operating temperatures: 4 – 40 °C
- Ends sealed with blanking caps, bags or aseptic connectors
- Assembled under ISO Class 8 conditions
- Assemblies to withstand gamma, X-ray or equivalent irradiation of 25 kGy (maximum 45kGy)
- Shelf life to be 12 – 24 months

Flow Path

SymphonX™



Competences

GEMÜ

GEMÜ

FUJIFILM

Disynth
biotechnologies

Quality management systems

- Manufactured in compliance with GMP
- Compliance with ISO 9001
- Manufacturing site registered with the FDA
- TSE/BSE-compliant manufacture
- Transport validation as per ASTM
- Cleanroom manufacturing and assembly in cleanroom class 8 according to EN ISO 14644-1

Product requirements

- ISO 10993 compatibility
- USP/EP-compliant purity
- Use of FDA-compliant plastics
- Radiation sterilization by X-ray



Design and production of individual products and solutions

Single-Use Diaphragm Valve

GEMÜ SU40 SUMONDO


Features

- Easy, application-friendly valve body assembly
- Impact-resistant and corrosion-resistant
- Integral optical position indicator
- Simplified installation thanks to the fastening clamp
- Optimized body replacement time



Single-Use Diaphragm Valve

GEMÜ SU40 SUMONDO | connection principle



Single-Use Diaphragm Valve

GEMÜ SU60 SUMONDO

Features

- Hermetic separation between medium and actuator
- Optimized assembly options
- Open/close function, positioner and process controller
- Force and speed are variably adjustable
- Optimized body replacement time
- Can be operated via eSy-Web or Modbus TCP web interface




maruX™

Continuous Monoclonal Antibody Manufacturing

Connected SymphonX™ case study

FUJIFILM Diosynth Biotechnologies



www.fujifilmdiosynth.com



Continuous mAb Production Facility


- 400 m² non-GMP suite (Billingham, UK)
- Functionally closed end-to-end disposable flow path
- Semi-continuous operation
 - Upstream: 500L Perfusion bioreactor and ATF cell retention device using Apollo™X perfusion cell line and FUJIFILM Irvine Scientific media
 - Downstream: intensified batch processing with SymphonX™
 - Intermediate break bags on load cells
 - Viral clearance and analytics are aligned with current fed-batch processes
 - Minimising 'Cost of Quality' and perceived Regulatory risk



Integrated USP and DSP generates > 15 kg purified mAb in 30 days

www.fujifilmdiosynth.com




Process Overview non-GMP



maruX™ Process Performance



- Highly productive perfusion process
 - Cell density 120×10^6 cells/mL $\pm 10\%$
 - Viability >90 %
 - Consistent product quality

- Robust semi-continuous downstream purification process
 - Expected clearance of residuals
 - Consistent product quality




maruX™ Buffer Management

- Point-of-use in-line dilution using SymphonX™
- Buffer Strategy
 - 5x buffer concentrates for buffers on all 7 unit operations
 - 4x buffer stock for high salt and formulation buffer
 - 1x buffers used for CIEX gradient
 - All pre-made and 0.2 mm filtered
- ~14,500 L liquid / week
 - 10,000 L Water
 - 1,600 L 5x buffer conc / week
 - 800 L 4x buffer conc / week
 - 2,100 L 1x buffer /week



www.fujifilmdiosynth.com

maruX™ Bioburden Control




www.fujifilmdiosynth.com

Drug Substance Quality Attribute		Target	Batch Process	Continuous Process
Concentration	(g/L)	35 – 45	37.2	39.5
SEC: Monomer peak	(%)	>87.5	97.3	98.4
Non Reduced CE-SDS	(%)	78.6 – 96.0	90.6	91.1
Reduced CE-SDS: Heavy Chain	(%)	58.5 – 71.5	65.0	64.2
Reduced CE-SDS: Light Chain	(%)	29.3 – 35.9	32.6	32.3
cIEF: pI		8.50 – 9.40	8.95	8.97
cIEF: Neutral peak	(%)	47.3 – 57.8	48.9	56.9
N-Glycan: GOF	(%)	65.8 – 80.4	73.1	74.8
Residual DNA	(pg/mg)	<100	2.09	<0.16
Residual HCP	(ng/mg)	<100	2.2	0.9
Residual Protein A	(ng/mg)	<20	Not determined	4.0

maruX™ Next step: Proposed GMP Operation

- **500 L perfusion bioreactor**
 - Constant feed into downstream ≥ 25 days
 - Residence time distribution mapping
- **Integrated downstream purification**
 - Flow path change out brackets lot
 - Traceability equivalent to current batch strategy
 - Viral clearance equivalent to current batch strategy



www.fujifilmdiosynth.com

Collaboration enabled advancement of technology platform

- Development of SymphonX™ through collaborative relationship with suppliers
 - Co-development of components
 - Delivery of system fit for purpose
- Successful demonstration of 500 L manufacturing scale continuous mAb production in a non-GMP facility
- Single-use mammalian GMP facility under construction in the UK for 2024 operations with 500 L GMP MaruX platform

maruX™

Apollo™ X
ADVANCED MAMMALIAN EXPRESSION™

SYMPHON X™

www.fujifilmdiosynth.com

Thank you for your attention

QUESTIONS?